3591 C-1

HIGHWAY RESEARCH REPORT

X-RAY DIFFRACTION



October, 1968

STATE OF CALIFORNIA

TRANSPORTATION AGENCY

DEPARTMENT OF PUBLIC WORKS

DIVISION OF HIGHWAYS

MATERIALS AND RESEARCH DEPARTMENT

RESEARCH REPORT

NO M&P 642002

DEPARTMENT OF PUBLIC WORKS

DIVISION OF HIGHWAYS

MATERIALS AND RESEARCH DEPARTMENT 5900 FOLSOM BLVD., SACRAMENTO 95819



October 1968

Final Report M&R No. 642992

Mr. J. A. Legarra State Highway Engineer

Dear Sir:

Submitted herewith is a research report titled:

X-RAY DIFFRACTION TECHNIQUES

TRAVIS W. SMITH Principal Investigator

M. L. McCAULEY and DR. J. L. POST Co-Investigators

Very truly yours,

JOHN L. BEATON

Materials and Research Engineer

ing section of the se

1 m

....

Side of the state of the state

REFERENCE: Post, J. L., "X-ray Diffraction Techniques", State of California, Department of Public Works, Division of Highways, Materials and Research Department, Research Report No. 642992, October 1968.

ABSTRACT: Various methods of quantitative analysis and identification of soil minerals by X-ray diffraction were investigated. The study disclosed that the needs of the Materials and Research Laboratory were best met by use of an oriented powder sample in a front pack to identify the non-clay minerals, and by use of strongly oriented coatings on glass slides for clay mineral identification. The slide samples were subjected to two heat treatments (325°C and 550°C), and to glycolation. Most of the clay minerals common in soils were identified by this method and numerous soil aggregates were quantified. The diffractograms of these samples were filed for reference.

The correlation of diffraction peak intensities to the amount of a mineral in an aggregate was studied by using powder samples in which the percent of the constituents was accurately known. A strong diffraction peak, or peaks, characteristic of the mineral was chosen, and the peak heights were plotted against the amount of mineral present. Calibration curves were then drawn from which the amount of mineral present can be determined.

The standard minerals used to establish calibration curves matched the condition of the minerals to be quantified. Changes in peak height caused by weathering were studied by making diffractograms of selected soil samples and observing the changes in reflection intensities. Experience was gained in the use of weathering factors when the standard minerals differed in degree of weathering from the soil minerals being quantified.

A detailed outline of the procedure for sample preparation, a list of auxiliary tests, and a table of $d(\mbox{\normalfont\AA})$ spacing for aggregate minerals were compiled.

The amount of iron contained in a sample was estimated by making diffractograms of chemically analyzed soil samples with variable iron content. The estimation was made by measuring the flourescence background intensity.

The flourescence background and absorption effects of iron and amorphous material in mineral samples were studied by making diffractograms with samples containing known amounts of iron and amorphous material. The calibration curves were nonlinear due to these absorption effects, but correction charts can be used.

KEY WORDS: Calibrations, clay minerals, minerals, quantitative analysis, samples, soils, X-ray diffraction.

tion of the second s The second secon

CONTENTS

		Page
Introducti	on	1
Summary		1
Sample Mou	nting Procedures	2
Packing Me	thods	4
Quantitati	ve Analysis Methods	5
Diffractio	n Peak Intensities	5
General Pr	ocedure	6
Auxiliary '	Tests	7
X-ray Diff	ractogram Analysis	7
Calibratio	n Curves	9
Diffractio	n Analysis Problems	15
References		20
	ILLUSTRATIONS	
Figure 1.	Effect of X-ray Tube Orientation and Resolution Slits on Diffraction Peak Intensity and Resolution	10
Figure 2.	Relationship Between Quartz Content and Four Diffraction Peak Intensities	11
Figure 3.	Relationship Between Mineral Contents, Muscovite, Hornblende, Serpentine, and Their Diffraction Peak Intensities	12
Figure 4.	Relationship Between Feldspar Content and 040 Diffraction Peak Intensities	14
Figure 5.	Fluorescent Background on Diffractogram Due to Iron Content of Aggregate Mineral Sample, 15° - 20° 20	16
Figure 6.	Effect of Hematite on Four Diffraction Peak Intensities, (110, 102, 112, 211)	17
Figure 7.	Effect of Amorphous Silic Acid, SiO2 on H2O, on Diffraction Peak Intensities	19

APPENDIX

						Page
Table 1.	d(Å)	Spacings	for	Aggregate Mine	erals	22-25
Calibrati	on Di	ffractogra	ams			26-31

The Make . The

orden otten.

Introduction

This study was made to evaluate the X-ray diffraction methods of mineral identification and quantitative analysis of soil and rock samples currently used by the Materials and Research Department. Also new procedures and techniques were to be evaluated, to assure that maximum benefit was being obtained from the equipment.

Dr. James L. Post, an Associate Professor of Civil Engineering at Sacramento State College, was retained as a consultant to accomplish the goals of this study. Dr. Post demonstrated for personnel of the Department proper preparation and treatment of samples, the proper use of the equipment and the interpretation of diffractograms. His extensive knowledge of the literature about X-ray diffraction and clay minerals provided many useful references.

As a result of this study, the Engineering Geology Unit is now using an entirely different X-ray diffraction procedure for the identification and quantitative analysis of the mineral composition of soils and rocks.

This research program was completed with the approval and financing of the State of California.

Summary

As a result of this research project, the following procedures have been adopted:

- 1. A strongly oriented sample is used. This sample is prepared by dense-packing the material using the front-pack method. Clay minerals are identified using two heat treatments (325°C and 550°C) and glycolation of a coating of the sample on a glass slide. This technique provides quick mineral identification and provides a fairly reliable quantitative analysis.
- Calibration curves showing particular diffraction peak intensities have been developed for use in determining the amount of a mineral present in a sample. The standard mineral used for making the calibration curves must be comparable in terms of crystallinity, chemical composition and degree of weathering to the unknown mineral in the sample. Uniformity in the fabrication of the front-pack samples is essential to obtain consistent results. Precise quantitative analysis of soil mineral aggregate will require auxiliary testing such as differential thermal analysis, microscopic examination or chemical analysis.

- 3. A systematic method for interpreting X-ray diffractograms was developed, and a list of the more common diffraction peaks was compiled and is included in this report.
- 4. The goniometer is now being run at 2° 2^{\bullet} per minute for routine scans and at 0.2° 2^{\bullet} if greater resolution is required.
- 5. To increase the resolution, the X-ray tube was turned from line focus to spot focus.
- 6. The chart recorder was adjusted so that one inch on the graph is equal to 5° 20. This permits visual indexing of the graph while the machine is still running.

In addition to developing the above described test procedure, certain other problems were studied. The effect of weathering on crystalline minerals in soil was particularly significant. The effects of iron fluorescence and amorphous minerals on the intensity of background radiation was studied. A method of estimating the amount of iron in a sample based on background intensity is discussed. Non-linearity of calibration curves, caused by X-ray absorption, is possible in certain samples and modification of the method of quantitative analysis based on linear relation-ships is necessary.

Sample Mounting Procedures

The identification of minerals which comprise rock and soil aggregates may be conducted by X-ray diffraction using two different types of sample preparations. Powder samples with varying degrees of random orientation obtained by one of several different packing procedures are used with the X-ray diffractometer, and a spindle or wedge containing powder samples are used with the X-ray powder camera. Oriented samples, consisting of powder packs or films of air-dried soil slurries on mounts, are used with the X-ray diffractometer and single-crystal mounts are used with the X-ray powder camera.

Quantitative analysis of mineral aggregations may be conducted by the use of oriented samples or samples with some degree of random particle orientation. The mineral analysis is ordinarily conducted in conjunction with supporting physical and chemical tests, such as clay-size content and equivalent iron sesquioxide content.

The degree of accuracy that may be obtained in the quantitative analysis of mineral aggregations by X-ray diffraction, with ordinary supporting tests, is largely dependent on the time available and the experience of the investigator. The mineral aggregation may be investigated as a whole or it may be separated into discrete particle-size ranges, such as the clay-size material. No one method should be adopted in exclusion of the others, rather the method should be used that best suits the degree of accuracy required and the time that is available for each case.

Rapid identification of aggregate minerals may be effected by using a front-pack powder sample and the mineral identification may be facilitated by firmly pressing the powder into the holder so that a strongly preferred orientation is obtained. Most soil minerals, other than quartz, tend to orient readily because of their cleavage characteristics, and because many soil aggregates are rather poorly crystalline the preferred orientation of the powder sample enhances the intensity of the diffraction peaks facilitating the identification and the quantitative analysis of the soil minerals.

The clay mineral suite is most readily identified by the process of preparing oriented coatings of the sample on glass slides and following the procedures given by Warshaw and Roy($\underline{1}$). The identification of small particles of hand-picked minerals from aggregates may usually be made by cementing a cluster of the particles on a glass slide and obtaining a diffraction pattern. Minerals on petrographic slides may be identified in the same manner according to Tatlock($\underline{2}$).

The <u>reliability</u> of quantitative analysis by X-ray diffraction is dependent on the consistency of sample preparation no matter which method of sample preparation is used. The <u>variation</u> in diffraction peak intensities of any mineral is largely dependent on the crystallite orientation and the bulk density of the powder sample, not to mention degree of crystallinity. Pressed pellets made from a powder sample give a preferred orientation to the sample comparable to a pressed front-pack powder sample, however, the sample bulk densities are probably more consistent. One disadvantage of the pressed pellet is the friability of pellets formed from framework silicates.

Random-orientation powder samples may be prepared both by front-pack and back-pack procedures, not to mention glass spindles for powder camera work. Comparable results may be obtained with each method, depending upon the experience and technique of the investigator. The powder sample may require a cover when using a General Electric X-ray diffractometer because the sample holder is placed in an upright position on the goniometer, whereas the Norelco sample holder is held in a horizontal position. A thin sheet of colloidon is sometimes used to cover the front-pack sample.

The back-pack procedure is often used because the method gives quite consistent results, even when used by different investigators. The procedure requires considerable time for the preparation of a randomly-oriented powder sample as does the randomly-oriented front-pack method, perhaps the larger part of an hour per sample in comparison to five minutes for the pressed front-pack method. The back-pack procedure does not give complete random orientation under many conditions of preparation, such as pressure on the back of the sample.

A list of advantages and disadvantages of different packing methods is included to give some idea of the many choices available.

Packing Methods

Front-Pack (pressed)

- Advantages: 1. Quick mineral identification and fairly reliable quantitative analyses.
 - 2. Moderate control of orientation (preferred orientation).
 - Intensification of necessary mineral identification peaks.
 - 4. Material passing a No. 200 sieve suffices.
- Disadvantages: 1. Some variation in degree of orientation and bulk density.
 - 2. Accuracy of results less than randomlyoriented front-pack.

Back-Pack

- Advantages: 1. Fairly consistent packing density.
 - 2. Consistent μ^* values determined by X-ray transmission.
- Disadvantages: 1. Irregular degree of orientation (feldspars strongly oriented).
 - 2. Major soil mineral identification peaks suppressed (orientation & lower density).
 - 3. Time-consuming preparation of samples.

Pellet

- Advantages: 1. Some control of preferred orientation and density.
 - 2. Quick preparation and mineral identification with fairly consistent peak heights.
- Disadvantages: 1. Pellets often friable causing diffraction trace distortion.
 - 2. Surface of material with high $\%\omega$ adheres to compactor.
 - 3. Sorting of mixture when preparing pellet.

Glass Slide (& porous filter method)

- Advantages: 1. Strong basal orientation for clay mineral identification.
 - 2. Amenable to quick heat treatment and fairly reliable results.
- Disadvantages: 1. Mineral sample thickness and consistency irregular.
 - 2. Only good for clay mineral determinations (water sorting).

Quantitative Analysis Methods

There are two ordinary methods of conducting quantitative analyses of mineral aggregations using an X-ray diffractometer: by determining the area under certain characteristic diffraction peaks, Klug and Alexander (3); by determining the heights of certain diffraction peaks, Tatlock (2); some combination of peak heights and contained areas, Schultz (4). In each case the intensity of diffraction is considered to be directly related to the amount of material present in the aggregation when the random-radiation background is accounted for.

The effect of physical weathering of minerals ordinarily has a far greater effect on diffraction peak intensities than any other factors, see page 422 ref. (3), except perhaps the degree of crystallinity of the aggregation when formed. Because of the number of minerals present in most soils and the physical weathering of the minerals, it is often necessary to depend upon diffraction peak heights for estimating the quantity of each mineral present. The large combination of diffraction peaks obscure one another, the peak heights are reduced by physical weathering, and the true radiation background is difficult to determine.

<u>Diffraction Peak Intensities</u>

The amount of a mineral in an aggregate may be correlated to diffraction peak intensities by several different methods. Ordinarily, a strong diffraction peak, or peaks, characteristic of the mineral is chosen which is not subject to undue interference from other mineral diffraction peaks and the peak height or enclosed area of the peak is related to the amount, by weight, of the mineral present. The calibration curves thus derived are not necessarily linear functions although many soil minerals appear to give linear correlations.

Standard minerals must be chosen which are representative of the minerals contained in the aggregate being investigated because the accuracy of an analysis depends directly on how well the standards are chosen, as indicated by Brindley, p. 501 (5).

The standard minerals may be derived from the aggregate samples or the minerals may be matched by materials from other sources. If other source material is used it <u>must</u> be comparable in terms of degree of crystallinity, chemical composition, and degree of weathering that has occurred.

The method of sample preparation is important because some soil minerals, such as clays and gypsum, are readily altered or destroyed by certain physical treatment such as over heating or over grinding. Powder samples consisting of a crystallite size of about 20μ tend to give best diffraction intensities and, perhaps, random orientation. It is not ordinarily necessary to grind material finer than 74μ (No. 200 sieve) for quantitative approximations if the pressed front-pack is used. A disc pulverizer may be used for an initial grind, followed by drygrinding of a split-out sample by mortar and pestle.

General Procedure

After the sample is received the following preliminary treatment should be done:

- 1. Air dry the sample and break down weakly cemented particles with rubber pestle.
- 2. (a) Coarse Aggregate Split out representative portion and reduce material to particles passing a No. 40 sieve using a disc pulverizer.
 - (b) Fine Aggregate (soils) Split out representative portion and separate the material passing the No. 10 sieve, weigh the material larger than 2 mm, and treat the larger material the same as the coarse aggregate if it comprises more than 10-15% of the total aggregate.
- 3. Determine the clay-size content of material passing the No. 10 sieve by the hydrometer analysis method. Clay-size material may be separated from the sample at this time as required. A cyclone separator may also be used.
- 4. Split out a portion, about 10-15 gms, of pulverized aggregate or soil and grind the material with a ceramic mortar and pestle. Reduce the aggregate and sandy or silty soil to material passing a No. 200 sieve (74) and reduce clayey soils to material passing a No. 140 sieve (0.1 mm). If clay-size material is separated from the sample it should be reduced to powder for mixing and handling purposes.
- 5. Put the reduced samples in sealed containers and mix thoroughly for further use.
- 6. Prepare samples of coarse-grind and soil aggregations for auxiliary physical-chemical tests (20-25 gms) where necessary.

Auxiliary Tests

It is seldom possible to conduct a very precise quantitative analysis of soil mineral aggregates without auxiliary tests. The tests required for different soils must be judged by the investigator. The tests may include:

- 1. Total carbonate content by gravimetric determination.
- 2. Total organic mineral content by heating or peroxide method.
- 3. Water-soluble content by leaching method.
- 4. Clay-size content by hydrometer analysis.
- 5. Ferric iron content by iodometric titration method or X-ray fluorescence background measurements.
- 6. Gypsum content by weight loss from dehydration (50° 100°C).
- Differential thermal analysis of clay minerals.
- 8. Heat treatments and "Glycolation" of oriented slides of clay minerals.
- 9. Optical investigation of coarse aggregates.
- 10. Specific chemical analyses such as Na+ cation content.

Less precise quantitative analyses of soils may be made by grinding a representative portion of the total air-dry sample with a disc pulverizer, splitting the coarse-grind material, and preparing a powder sample in the same manner as with soil fines. In any case, the moisture content of the air-dry soil sample (-#10 size) is required, the oven-dry temperature not to exceed 50°C.

X-ray Diffractogram Analysis

It is convenient to make the diffractogram so that one inch on the chart paper is equal to $5^{\circ}2\theta$ because sufficiently precise 2θ values may then be derived by visual estimate. The goniometer probably should not be run faster than $2^{\circ}2\theta$ per minute because the recorder does not have a sufficiently rapid response to keep up with a rapidly varying signal input. Tatlock (2) considers a $2^{\circ}2\theta$, drive rate to give results nearly comparable to $1^{\circ}2\theta$, which is more commonly used for precise mineral analysis. Use machine recorder constants of 500 CPS, TC = 2, with the G.E. XRD-5 unit.

The diffractogram should be precisely labelled while it is being made and the diffraction 20 angles and matching "d" spacings may be marked on the chart at that time if convenient. For a complete mineral identification all real peaks above the radiation background should be marked. For quantitative analysis only certain diffraction peaks will be used. A list of the more common diffraction peaks which indicate the presence of different aggregate minerals is given in Table 1, the peaks commonly used for quantitative analysis being marked by an asterisk.

The identification of the aggregate minerals is done in the following manner (3° - 63° 2θ):

- 1. The quartz peaks are marked. Often 6 of the 10 more prominent peaks are visible. When less than 10 percent quartz is present only 2 peaks may show, 4.26Å and 3.343Å. The 3.343Å quartz peak will often appear when as little as 3 percent quartz is present.
- 2. The feldspar peaks are marked. When more than about 10 percent feldspar is present it is necessary to make a second scan from 26.0° 29.0°20 at 2000 CPS, TC = 1, with a goniometer drive rate of 0.2° or 0.4°20 per minute, to identify the predominant feldspars so that the feldspar peaks may be properly identified. The 040 diffraction peaks are the main sources of feldspar identification.
- The amphibole and pyroxene peaks are marked. The most common amphibole is hornblende with primary identification peaks at 8.4Å and 3.12Å. Several of the amphiboles, constituting a solid state series, have primary peaks near 8.4Å 'd' spacing. The most common pyroxene is augite with primary identification peaks near 2.99Å, 2.94Å, and 1.62Å. Augite chemical composition is also variable.
- 4. The carbonate and sulfate peaks are marked. The primary identification peaks for the carbonates are: Calcite 3.03Å; Dolomite, 2.89Å; and Magnesite, 2.74Å. Aragonite is metastable in soils and not too common. The only common sulfate, gypsum, has a primary identification peak at 7.56Å, but may exist as hemihydrate in dessicated soil. Anhydrite is not common in soils nor is 7 CaSO4, the latter being found in gypsiferous soils oven-dried at 110°C, Deer, et al, vol. 5, p. 207 (6).
- 5. The mica, chlorite, and serpentine peaks are marked. The main diffraction peak for micas occurs at 10.0Å, chlorites at 7.0-7.2Å, and serpentines at 7.36Å. In soils these minerals may be identified along with the clay mineral suite.
- 6. The less common non-clay mineral peaks are marked. The main diffraction peaks which suggest the presence of these minerals are given in Table 1. If any of these minerals are believed to be present the entire diffraction pattern listing should be checked for those minerals. Halite, magnetite, goethite, zeolites, and cristobalite often occur as minor soil components (< 5%).
- 7. The clay mineral peaks are marked. The known clay mineral peaks on the powder sample diffractogram are marked and the oriented slide technique is used generally according to the procedure of Warshaw and Roy (1) to identify the clay minerals.

The relative intensities given for X-ray diffraction patterns for different minerals by ASTM Powder Diffraction File and other sources of information are supposed to be for randomly oriented powder samples. The intensity values may be used as a rough guide in mineral identification but diffraction peak intensities obtained in the laboratory are commonly quite different, depending upon the cleavage habits of the minerals and the degree of preferred orientation achieved. Of course, the crystallite size and degree of crystallinity may cause some variations in peak intensities.

Ordinarily powder sample scans are made from 3° to 63° 29 and oriented slide scans are made from 3° to 21° 20 at 500 CPS, TC = 2. If it is necessary to identify very poorly crystalline clays a count of 200 CPS may be used. Medium resolution Soller slits are normally used but the high-resolution slits with the 0.05° receiver slit, better resolve the feldspar peaks from 26° to 29° 20 at 2000 CPS, TC = 1 using a goniometer drive rate of 0.2° 20 per minute. Much better peak resolution may be obtained by turning the X-ray tube to center port for "spot" focus and using a 0.05° receiver slit, as shown in Figure 1.

Calibration Curves

Calibration curves may be made for crystalline aggregate minerals by relating selected diffraction-peak intensities to the amount of each mineral present in the aggregate. Calibration curves are derived by mixing together known amounts of soil minerals, or by using previously analyzed mixtures, and making diffractograms of the mineral calibration peaks. The peak heights are then plotted against the amount of mineral present. In some cases the areas under the diffraction peaks are used for quantitative analysis of different minerals. Diffraction peaks commonly used to make calibration curves are marked with asterisks in Table 1.

Quartz is the only major crystalline constituent of rock and soil aggregates that has very poor cleavage, thus, the powder crystallite orientation is quite variable when using a pressed powder sample. By using four calibration peaks for quartz, as shown in Figure 2, and taking the average value a close approximation of quartz content (by weight) may be obtained. The four calibration curves may be considered linear and represent approximately a random powder orientation. Curves representing relatively unweathered quartz may be derived by using diffraction intensity data obtained with a permaquartz slab, the peak heights being plotted as 100% quartz as shown in Figure 2.

Calibration curves for three common aggregate minerals, muscovite, hornblende, and serpentine, are given in Figure 3. The powder samples which were prepared consisted of relatively unweathered material, the muscovite occurring as phyllite, the hornblende as finely crystalline rock, and the serpentine as antigorite. All three of the minerals orient readily and thus




Figure 1. Effect of x-ray tube orientation and resolution slits on diffraction peak intensity and resolution.

(a) outer peaks: line source, 200 CPS; inner peaks: spot source, 500 CPS. 0.05° receiver slit. (b) outer peaks: 0.2° slit, middle peaks: 0.1° slit, inner peak: 0.05° receiver slit. Spot source.




Figure 2. Relationship between quartz content and four diffraction peak intensities, using permaquartz. (G. E. SPG-3 receiver, Cu/Ni, 500 CPS, TC 2, slits: 10-.20 mr)




Figure 3. Relationship between mineral contents, muscovite, hornblende, serpentine, and their diffraction peak intensities.

(G. E. SPG-3 receiver, Cu/Ni, 500 CPS, TC 2, slits: 10-.20 mr)

give quite reliable results when using the pressed-powder packing technique.

The feldspars usually occur as perthite, a sodic plagioclase in microcline, and require special treatment for detailed identification. Feldspars possess perfect cleavage and the powder crystallites tend to orient so that the 0k0 diffraction peaks are enhanced. Calibration curves may be made using feldspar 040 diffraction-peak intensity measurements. Because of the intensity of the diffraction peaks and the close spacings of the feldspar 040 peaks, it is necessary to scan at a slow rate $(0.2^{\circ}~2\theta~\text{minute})$ and to change count rate from 500 CPS to 2000 CPS when more than about 10 percent feldspar is present. Figure 4 shows a typical calibration curve for total feldspar content $(\Sigma~040~\text{peak~heights})$.

The standard minerals which are used to establish calibration curves must match the condition of the minerals which are to be quantified, especially with respect to degree of crystallinity and chemical composition. If the minerals in the aggregate only vary from the standard minerals in the degree of physical weathering, the most common situation, a weathering factor must be introduced or additional standard minerals must be secured from the aggregate or its equivalent. When many rock and soil aggregates with varying degrees of physical weathering are to be considered it is necessary to establish weathering factors to account for the decline in intensity of the mineral diffraction peaks.

Soil and rock minerals vary in their resistance to physical and chemical weathering so that each mineral will have a distinct weathering range. Some minerals are very stable and show little effect from weathering. Quartz, for example, shows a peak diffraction intensity for dune sand or extremely weathered tropical soil that is about two-thirds that of unweathered material. Chalcedony (cryptocrystalline a-quartz) shows about the same peak intensity as extremely weathered quartz. A maximum weathering factor (multiplier) of 1.50 might be used for quartz to adjust to standard calibration curves for unweathered quartz, and, of course, there are also degrees of weathering to consider.

Most soil and rock aggregates are crystalline, excepting unweathered igneous extrusives which contain glassy material, and because of this crystallinity the sum of the minerals, represented by diffraction peaks, must total nearly 100 percent when including organic material and water content. Knowing the rate and degree of weathering possible for different minerals a reasonable estimate of the aggregate mineral content may be made. Some soil minerals are very subject to chemical weathering and do not appear to require weathering factors, including amphiboles, pyroxenes, and serpentines. Either they are present as formed or they are chemically altered. Calibration curves for authigenic soil minerals, such as sedimentary chlorite, must be derived partly from allied tests. As an example,




Figure 4. Relationship between feldspar content and 040 diffraction peak intensities. (G. E. SPG-3 receiver, Cu/Ni, 2000 CPS, TC 1, slits: 10-.10 mr)

certain extremely weathered tropical soils were found to consist of only quartz and kaolinite with a trace of goethite, Post (7). The hydrometer analysis established the clay-size content (kaolinite) and the equivalent ferric iron content was found chemically (giving goethite), the remaining mineral then being quartz. By deliberately choosing simple two, three, or four-mineral systems the weathering stage of aggregates of an area may thus be determined.

Some comparison of the effect of weathering on different minerals may be shown by considering perthite in relation to quartz. Where extremely weathered quartz may require a weathering factor of 1.5, extremely weathered perthite in rock may require a factor of 2 and in soil a factor of 3. Whereas micas are very stable minerals in aggregates with little change in crystallinity from physical weathering, clay minerals such as kaolinite and montmorillonite may appear to be nearly amorphous after extreme physical weathering. Illite (hydromica) is ordinarily considered to be formed by the chemical weathering of micas, part of the K⁺ ions being removed from the interlayer position. The K⁺ ions are weakly bonded, being held in 12-fold coordination by the silica Tetrahedra. Illite and kaolinite in soil aggregates ordinarily show about the same diffraction intensities for first order (001) basal spacings.

Diffraction Analysis Problems

Under certain conditions the total iron in soil and rock aggregates may be estimated by measuring the fluorescence background intensity on the diffractogram, Figure 5, and when the iron is known to be contained in one mineral the amount of mineral contained may also be estimated, as shown by Tatlock (2). Several common aggregate minerals contain variable amounts of iron, including hornblende, augite, and biotite, and many soils also include clay minerals which contain iron. The common soil also include clay minerals which contain iron. The common somineral goethite, Fe203.H20, is not very susceptible to X-ray diffraction and ordinarily must be quantified by chemical The weakness of the X-ray reflections are presumably due to fine grain size rather than poor structural order, Weaver, et al (8). It may be seen, Figure 5, that the background caused by goethite appears to include not only radiation scatter caused by iron fluorescence but also scatter caused by the microcrystalline nature of goethite, and possibly by the presence of amorphous material.

Iron is the only relatively abundant common rock-forming element whose fluorescence under copper radiation effects appreciable increases in background radiation. Not only does the fluorescence background caused by iron tend to cover other mineral diffraction peaks, Figure 6, but because iron is a strong radiation absorber ($\mu = 324$), peaks of the weakly absorbing component appear weaker and those of the strongly absorbing component appear stronger, Klug and Alexander (3).




Figure 5. Fluorescent background on diffractogram due to iron content of aggregate mineral sample, 15° - 20° 20. (Gas flow counter, 500 CPS, TC 2)




Figure 6. Effect of hematite on four diffraction peak intensities, sum of quartz 110, 102, 112, and 211 reflections. (Gas flow counter. 500 CPS, TC 2)

The effects of iron in the X-ray diffraction analysis of aggregates may be adjusted by the use of correction charts as Tatlock (2) has proposed or the effects of iron may be negated by using an iron or chromium X-ray tube target. Unfortunately, the iron and chromium X-ray tubes do not give as useful diffraction patterns as the copper tube.

Amorphous minerals raise the diffraction background intensity in varying amounts, probably in some relationship to the mass absorption coefficient of the amorphous material present. Silicic acid, SiO2·nH2O, gives a very low background intensity (μ * = 29.7) whereas micro-crystalline goethite gives a very high intensity (μ * = 208.5).

In comparing diffraction patterns of aggregate powders, absorption effects are present that cause calibration curves to be non-linear. Whenever mass absorption coefficients of any two predominant aggregate minerals vary considerably the calibration curves may appear as shown in Figure 7. This does not hinder the use of the analysis procedure but it must be remembered that the linear relationship given by Klug and Alexander (3) may not be used without some modification.




Figure 7. Effect of amorphous silicic acid, SiO₂*nH₂O, on diffraction peak intensities of gypsum. (Gas flow counter, 500 CPS, TC 2)

REFERENCES

- Warshaw, C. M., and Roy, R. "Classification and a Scheme for the Identification of Layer Silicates," Geol. Soc. of America <u>Bulletin</u>, Vol. 72, No. 10, Oct. 1961, p. 1455-1492.
- Tatlock, D. B. "Rapid Modal Analysis of Some Felsic Rocks from Calibrated X-ray Diffraction Patterns," Geol. Survey <u>Bulletin</u> 1209, Wash., D.C., 1966.
- 3. Klug, H. P., and Alexander, L. E. X-ray Diffraction Procedures for Polycrystalline and Amorphous Materials, John Wiley and Sons, Inc., N.Y., 1954, p. 410-438.
- 4. Schultz, L. G. "Quantitative Interpretation of Mineralogical Composition from X-ray and Chemical Data for the Pierre Shale," Geological Survey Prof. Paper 391-C., Wash., D.C., 1964.
- 5. Brown, G. (ed) The X-ray Identification and Crystal Structures of Clay Minerals, Mineralogical Society, London, 1961, p. 489-514.
- 6. Deer, W. A., Howie, R. A., and Zussman, J. Rock-Forming Minerals, Vols. 1-5, John Wiley and Sons, Inc., N.Y., 1962.
- 7. Post, J. L. Engineering Characteristics of Some Southeast

 Asian Soils, Technical Report No. AFWL-TR-67-143,

 Air Force Weapons Laboratory, Kirtland AFB, New Mexico,
 May 1968.
- 8. Weaver, C. E., Wampler, J. M., and Pecuil, T. E. "Mossbauer Analysis of Iron in Clay Minerals," <u>Science</u>, Vol. 156, April 1967, p. 504-508.

APPENDIX

 $d(\mathring{A})$ Spacings for Aggregate Minerals Calibration Diffractograms

Table 1

d (Å) Spacings for Aggregate Minerals


	<u>d(Å)</u>	<u>I</u>	Mineral	<u>Notes</u>
*	15.0-12.5	vs	Mont. Group (Smectite)	varying degrees of hydration
	14.5-13.5	vs	Chlorite	
*	14.5-13.5	vs	Vermiculite	fully hydrated
*	12.5-12.0	vs	Vermiculite	partially hydrated
*	10.5	vs	Attapulgite	(palygorskite)
*	10.1	S	Illite	
*	10.1-9.9	vvs	Micas	muscovite, biotite, phlogopite
*	9.3	s ·	Talc	
*	8.5-8.4	vs	Amphibole	hornblende, actinolite
*	7.56	vs	Gypsum	
ж̂	7.36	vs	Serpentine	antigorite, chrysotile
*	7.5-7.2	wm	Halloysite	metahalloysite (2H ₂ O)
*	7.15	s	Kaolinite	
*	7.2-7.0	vs	Chlorite	iron-rich chlorites to 6.78Å
	6.65	m	Orthoclase	"high" K-feldspar
	6.46	m	Microcline	"low" K-feldspar
	6.55-6.38	m	Plagioclase	anorthite to albite
*	5.61-5.57	mw	Analcite	
*	4.85	S	Gibbsite	
,	4.79-4.69	ms	Chlorite	
*	4.50-4.45	m	Clay 020-040	mainly kaolinite & smectites
*	4.37	m	Gibbsite	
	4.27	s	Gypsum	
	4.26	s	Quartz	and K-feldspar 4.24-4.21Å


^{(*} Peaks used for quantitative analysis)


	<u>d(Å)</u>	I	<u>Mineral</u>	Notes
*	4.18	vs	Goethite	
*	4.04	vs	Cristobalite	and plagioclase peak
	4.04-4.02	mW	Plagioclase	albite to anorthite
	3.97-3.84	mw	Feldspars	
	3.83-3.79	s	K-feldspar	microcline & orthoclase
	3.77-3.74	ms	Plagioclase	
	3.66-3.61	ms	Plagioclase	
	3.66-3.64	s	Serpentine	
	3.57	ms	Kaolinite	
	3.59-3.52	s	Chlorite	
	3.48-3.46	m	Feldspars	
	3.43	S	Analcite	
	3.37-3.36	vvs	Micas	muscovite and biotite
	3.343	vvs	Quartz	
	3.33	vs	Orthoclase	
*	3.245-3.225	vvs	Microcline	and orthoclase
*	3.21-3.19	vvs	Plagioclase	albite to anorthite
	3.18	vvs	Plagioclase	only albite and oligoclase
	3.17-3.16	s	Plagioclase	only andesine to anorthite
	3.13-3.11	wm	Plagioclase	only andesine to anorthite
	3.13-3.09	S	Amphibole	hornblende and actinolite
	3.10	vs	Talc	
	3.09	vs	Jarosite	
	3.06	s	Gypsum	. 26
*	3.03-3.02	vvs	Calcite	No. He
•	3.01	vs	Alunite	张文


. "	<u>d(Å)</u>	I	<u>Mineral</u>	<u>Notes</u>
	3.01-2.95	ms	Feldspars	131 Peak
	2.99	s	Augite	
	2.94-2.91	mW	Feldspars	
*	2.89	vs	Dolomite	
	3.87-2.82	mw	Chlorite	
	2.86-2.82	mw	Feldspars	131 peak
	2.82	vs	Halite	
-	2.75	S	Ilmenite	
	2.74	vs	Magnesite	
	2.71-2.69	s	Amphibole	hornblende & actinolite
	2.69	sb	Hematite	and goethite
.*	2.56	S	Augite	
	2.56-2.51	m	Feldspars	and micas
	2.53	s	Magnetite	
	2.51	s	Augite	and hematite
	2.50-2.44	шМ	Feldspars	albite to labradorite
. '	2.48	m .	Cristobalite	
*	2.458	mw	Quartz	
*	2.33	m	Kaolinite	
*	2.282	m₩	Quartz	and calcite
٠.	2.19	ms	Dolomite	
	2.10	m₩	Calcite	
	2.01-1.99	m	Micas	muscovite & biotite
	1.99	s	Halite	
	1.91	шW	Calcite	
	1.87	mw	Calcite	


	<u>d(Å)</u>	I	<u>Mineral</u>	<u>Notes</u>
*	1.817	mW	Quartz	and dolomite
	1 482-1.80	mw	Feldspars	
	1.78	dwm	Dolomite	
	1.62	mw	Magnetite	
	1.56	ms	Septichlorites	esp. Chamosite (060)
	1.55-1.53	W	Chlorites	060 peak
*	1.541	mw	Quartz	211 peak
	1.54	mw	Biotite	060 peak
	1.54-1.53	W	Vermiculite	and trioct. smectite (060)
	1.53	wm	Serpentine	and talc (060)
	1.53-1.52	W	Phlogopite	060 peak
	1.52-1.49	W	Smectite	dioct. 060 peak
	1.50	m	Muscovite	and illite & attapulgite (060
	1.49	m	Kaolinite	and halloysite (060)

















1-55-3	8 8	
<u> </u>		
	10.0%	
	10.0	
	N N	
	3.E	70
	Z.V.	17
	-HORNBLENDE	. Ug.
	SERPENTINE 10.0%	MUSCOVITE 2,5%
	10,0	2
20	ш	38
40.0	TIN	₹
	E	
VAS	ERPE	
d	γ ∷	
1 译.		
+	1 /	3
	/ 8	
	\\$	
1		
	HH CA	LIBRATION T
		<u>i</u>

MUSCONITE 15:097	- HARM BLENDE 15:0%	Musequite 156%		45% FELDSPAR 10% QUARTZ 15% AUSCOVITE
<	SERPENTINE 150 %		10.1 % FELDSPAR 74-9 % QUARTZ	±an±1.
	- FELSSPAK 4567		FNDE.	TO STATE OF THE ENDE
	18/2		NATARALIAN 121	5 # CAL BRA HOL
	5.71 CALLE	RATION	TOO CPS TC=2 0.2° REC. SL. R° 26/MIN.	500 CPS TC=2 OZEREC, SL.

Alexandra de deservadores.